




RELB (Ab-573) Antibody

Product Code	CSB-PA697302
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01201
Immunogen	Peptide sequence around aa. 550~554 (L-L-S-P-G) derived from Human RelB.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total RelB protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	I-Rel
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RELB

Image



Western blot analysis of extracts from HUVEC, 293 and HeLa cells using RelB(Ab-573) Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.