




GRIA1 (Ab-863) Antibody

Product Code	CSB-PA226902
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42261
Immunogen	Peptide sequence around aa.861~865 (R-N-S-G-A) derived from Human GluR1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total GluR1 protein.
Tested Applications	ELISA,IF;IF:1:100-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	GLR1; GLUH1; GRIA1; GluR-1; GluR-A
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GRIA1

Image



Immunofluorescence staining of methanol-fixed HeLa cells using GluR1(Ab-863) Antibody.