




MARCKS (Ab-162) Antibody

Product Code	CSB-PA266770
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29966
Immunogen	Peptide sequence around aa. 160~164 (K-K-S-F-K) derived from Human MARCKS.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total MARCKS protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	MACS; MARCS; PKCSL; PRKCSL; Protein kinase C substrate
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MARCKS

Image



Western blot analysis of extracts from HeLa and HT29 cells using MARCKS(Ab-162) Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.