





CTTN (Ab-466) Antibody

Product Code	CSB-PA095187
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14247
Immunogen	Peptide sequence around aa.464~468 (P-V-Y-E-T) derived from Human CORTACTIN.
Raised In	Rabbit
Species Reactivity	Mouse
Specificity	The antibody detects endogenous level of total Cortactin protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:1000, IHC: 1:50-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	Amplaxin; CTTN; EMS1; SRC8;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CTTN

Image



Western blot analysis of extracts from 3T3 and C2C12 cells using Cortactin(Ab-466) Antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using Cortactin(Ab-466) Antibody(left) or the same antibody preincubated with blocking peptide(right).