




PPP1R14A (Ab-38) Antibody

Product Code	CSB-PA794704
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96A00
Immunogen	Peptide sequence around aa.36~40 (R-V-T-V-K) derived from Human CPI17a.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total CPI17a protein.
Tested Applications	ELISA,IF;IF:1:100-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	CPI; CPI-17; CPI17; PP14A;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPP1R14A

Image



Immunofluorescence staining of methanol-fixed HeLa cells using CPI17a(Ab-38) Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.