





NFKBIE (Ab-22) Antibody

Product Code	CSB-PA148874
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00221
Immunogen	Peptide sequence around aa. 20~24 (I-E-S-L-R) derived from Human IκB-e.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous level of total IκB-e protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:1000, IHC: 1:50-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	I-kappa-B-epsilon
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NFKBIE

Image



Western blot analysis of extracts from HeLa, 293 and HepG2 cells using IκB-e(Ab-22) Antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using IκB-e(Ab-22) Antibody(left) or the same antibody preincubated with blocking peptide(right).

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.