




OPRM1 (Ab-375) Antibody

Product Code	CSB-PA148947
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35372
Immunogen	Peptide sequence around aa.373~377 (H-P-S-T-A) derived from Human Opioid Receptor.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total Opioid Receptor protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	MOR; MUOR1; Mu-type opioid receptor; OPRM; OPRM1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	OPRM1

Image



Western blot analysis of extracts from Rat brain and Mouse brain tissue using Opioid Receptor(Ab-375) Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.