




PROX1 Antibody

Product Code	CSB-PA298119
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92786
Immunogen	Peptide sequence around aa.723~727 (E-I-F-K-S) derived from Human Prospero homeobox protein 1.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total Prospero homeobox protein. 1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	PROX1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PROX1

Image



Western blot analysis of extracts from A431 cells using Prospero homeobox protein 1 Antibody and the same antibody preincubated with blocking peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.