





FLT4 Antibody

Product Code	CSB-PA277523
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35916
Immunogen	Peptide sequence around aa.1279~1283 (L-A-S-E-E) derived from Human VEGFR-3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total VEGFR-3 protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:200,IF:1:100-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	PCL; FLT4; FLT41; LMPH1A;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FLT4

Image



Western blot analysis of extracts from HUVEC cells using VEGFR-3 Antibody (right) and the same antibody preincubated with blocking peptide (left).



Western blot analysis of extracts from Mouse liver tissue using VEGFR-3 Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.