



ORF4 Antibody

Product Code	CSB-PA133047
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15184
Immunogen	Peptide sequence around aa.9~13 (D-A-R-E-Q) derived from Cucumber necrosis virus RNA silencing suppressor.
Raised In	Rabbit
Species Reactivity	Virus
Specificity	The antibody detects endogenous level of total RNA silencing suppressor protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Cucumber necrosis virus (CNV)
Target Names	ORF4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.