



Fbxl4 Antibody

Product Code	CSB-PA064959
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9VY46
Immunogen	Peptide sequence around aa.1~5 (M-S-L-L-A) derived from Fruit fly F box and leucine-rich-repeat gene 4.
Raised In	Rabbit
Species Reactivity	Drosophila melanogaster
Specificity	The antibody detects endogenous level of total F box and leucine-rich-repeat gene 4 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	FBL4; FBL5; FBXL4
Product Type	Polyclonal Antibody
Immunogen Species	Drosophila melanogaster
Target Names	Fbxl4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.