




# MTOR (Ab-2481) Antibody

<b>Product Code</b>	CSB-PA054657
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P42345
<b>Immunogen</b>	Peptide sequence around aa.2479~2483 (I-H-S-F-I) derived from Human mTOR.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Specificity</b>	The antibody detects endogenous level of total mTOR protein.
<b>Tested Applications</b>	ELISA,IHC;IHC:1:50-1:100
<b>Form</b>	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	FRAP; FRAP1; FRAP2; RAFT1; Rapamycin target protein
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	MTOR

## Image



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using mTOR(Ab-2481) Antibody(left) or the same antibody preincubated with blocking peptide(right).

**Usage** For Research Use Only. Not for use in diagnostic or therapeutic procedures.