





CCNB1 (Ab-147) Antibody

Product Code	CSB-PA973969
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14635
Immunogen	Peptide sequence around aa.145~149 (A-F-S-D-V) derived from Human Cyclin B1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total Cyclin B1 protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:1000,IHC:1:50-1:200,IF:1:100-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CCNB1


Image



Western blot analysis of extracts from HepG2 and HT29 cells using Cyclin B1(Ab-147) Antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using Cyclin B1(Ab-147) Antibody(left) or the same antibody preincubated with blocking peptide(right).



Immunofluorescence staining of methanol-fixed HeLa cells using Cyclin B1(Ab-147) Antibody.