




BCL2L1 (Ab-62) Antibody

Product Code	CSB-PA923588
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07817
Immunogen	Peptide sequence around aa.60~64 (A-D-S-P-A) derived from Human BCL-XL .
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total BCL-XL protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	Apoptosis regulator Bcl-X
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BCL2L1

Image



Western blot analysis of extracts from 293 and HeLa cells using BCL-XL(Ab-62) antibody.