




# POU5F1 Antibody

<b>Product Code</b>	CSB-PA694841
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q01860
<b>Immunogen</b>	Peptide sequence around aa.232~236 (R-K-R-T-S) derived from OCT-4
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Specificity</b>	The antibody detects endogenous level of total OCT-4 protein.
<b>Tested Applications</b>	ELISA,WB;WB:1:500-1:1000
<b>Form</b>	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Otc3; OTF3; OTF4; OTF-3;
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	POU5F1

**Image**


Western blot analysis of extracts from human Umbilical cord mesenchymal stem cell using OCT-4 Antibody.

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.