





PRKCZ (Ab-410) Antibody

Product Code	CSB-PA927558
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05513
Immunogen	Peptide sequence around aa. 408~412 (T-S-T-F-C) derived from PKCz
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total PKCz protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:1000,IHC:1:50-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	PKC2; PKC-ZETA; PRKCZ
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRKCZ

Image



Western blot analysis of extract from 293 cells and HeLa using PKCz(Ab-410) Antibody.



Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue using PKCz(Ab-410) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.