




PTK2 (Ab-397) Antibody

Product Code	CSB-PA910905
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05397
Immunogen	Peptide sequence around aa. 395~399 (D-D-Y-A-E) derived from Human FAK.
Raised In	Rabbit
Species Reactivity	Human,Rat
Specificity	The antibody detects endogenous levels of total FAK protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	FADK; PTK2; FAK1; FRNK; PTK2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PTK2

Image



Western blot analysis of extract from HepG2, 293, C6, HUVEC and HT29 cells using FAK(Ab-397) Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.