




GRIA1 Antibody

Product Code	CSB-PA199388
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42261
Immunogen	Peptide sequence around aa. 831~835 (K-S-R-S-E) derived from Human GluR1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total GluR1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	GLR1; GLUH1; GRIA1; GluR-1; GluR-A
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GRIA1

Image



Western blot analysis of extracts from mouse brain and using GluR1 Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.