




GRIA1 Antibody

Product Code	CSB-PA199388
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42261
Immunogen	Peptide sequence around aa. 831~835 (K-S-R-S-E) derived from Human GluR1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total GluR1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	GLR1; GLUH1; GRIA1; GluR-1; GluR-A
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GRIA1

Image



Western blot analysis of extracts from mouse brain and using GluR1 Antibody.