




TFE3 Antibody

Product Code	CSB-PA104708
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19532
Immunogen	Peptide sequence around aa.303~307 (L-D-V-Y-S) derived from Human TFE3.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total TFE3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	BHLHE33; TFEA; RCCP2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TFE3

Image



Western blot analysis of extracts from 293 cells using TFE3 Antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.