




FLI1 Antibody

Product Code	CSB-PA282905
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01543
Immunogen	Peptide sequence around aa.251~255(S-K-N-T-E)derived from Human FLI1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total FLI1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	EWSR2; SIC-1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FLI1

Image



Western blot analysis of extracts from 293 and HeLa cells using FLI1 Antibody.