




NCAM1 Antibody

Product Code	CSB-PA960581
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13591
Immunogen	Peptide sequence around aa.850~854(Q-T-K-E-N)derived from Human CD56(NCAM).
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total NCAM protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	MSK39; NCAM1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NCAM1

Image



Western blot analysis of extract from rat brain, rat hippocampus and mouse brain using CD56(NCAM) Antibody.