




HNRNPA1 (Ab-192) Antibody

Product Code	CSB-PA737011
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09651
Immunogen	Peptide sequence around aa.190~194(S-S-S-Q-R)derived from Human hnRNP A1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total hnRNP A1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	HNRPA1; HNRPA1L3; hnRNP A1; hnRNP-A1; MGC102835
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HNRNPA1

Image



Western blot analysis of extract from 293 cells using hnRNP A1(Ab-192) Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.