




HNRNPA1 (Ab-192) Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA737011 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P09651 |
| Immunogen | Peptide sequence around aa.190~194(S-S-S-Q-R)derived from Human hnRNP A1. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Specificity | The antibody detects endogenous levels of total hnRNP A1 protein. |
| Tested Applications | ELISA,WB;WB:1:500-1:1000 |
| Form | Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Purification Method | Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide. |
| Clonality | Polyclonal |
| Alias | HNRPA1; HNRPA1L3; hnRNP A1; hnRNP-A1; MGC102835 |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | HNRNPA1 |

Image



Western blot analysis of extract from 293 cells using hnRNP A1(Ab-192) Antibody.