




TH (Ab-31) Antibody

Product Code	CSB-PA695636
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07101
Immunogen	Peptide sequence around aa.29~33(V-T-S-P-R)derived from Mouse Tyrosine Hydroxylase.
Raised In	Rabbit
Species Reactivity	Mouse,Rat
Specificity	The antibody detects endogenous levels of total Tyrosine Hydroxylase protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	TYH; DYT14; DYT5b; TH;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TH

Image



Western blot analysis of extract from rat brain and mouse brain using Tyrosine Hydroxylase(Ab-31) Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.