




Gripap1 Antibody

Product Code	CSB-PA210251
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9JHZ4
Immunogen	Peptide sequence around aa. 814-818(Q-E-I-V-R) derived from Rat GRASP.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total GRASP protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide
Clonality	Polyclonal
Alias	Gripap; Grasp; GRIP1-associated protein
Product Type	Polyclonal Antibody
Immunogen Species	Rattus norvegicus (Rat)
Target Names	Gripap1

Image



Western blot analysis of extracts from rat brain and mouse brain tissue using GRASP Antibody