




Gripap1 Antibody

Product Code	CSB-PA210251
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9JHZ4
Immunogen	Peptide sequence around aa. 814-818(Q-E-I-V-R) derived from Rat GRASP.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total GRASP protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide
Clonality	Polyclonal
Alias	Gripap; Grasp; GRIP1-associated protein
Product Type	Polyclonal Antibody
Immunogen Species	Rattus norvegicus (Rat)
Target Names	Gripap1

Image



Western blot analysis of extracts from rat brain and mouse brain tissue using GRASP Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.