




CD34 Antibody

Product Code	CSB-PA269036
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28906
Immunogen	Peptide sequence around aa.377~381(Q-H-V-V-A) derived from Human CD34.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total CD34 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	CD34
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CD34

Image



Western blot analysis of extract from HUVEC cells and Hematopoietic Stem Cells from Umbilical Cord Blood using CD34 Antibody