





CD44 Antibody

Product Code	CSB-PA644696
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16070
Immunogen	Peptide sequence around aa.734~738(Q-N-V-D-M) derived from Human CD44.
Raised In	Rabbit
Species Reactivity	Human,Rat
Specificity	The antibody detects endogenous level of total CD44 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:1000,IF:1:100-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	LHR; MDU2; MDU3; MIC4;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CD44

Image



Western blot analysis of extract from HeLa cells using CD44 Antibody



Immunofluorescence staining of methanol-fixed NCI-h446 cells showing cytoplasmic staining using CD44 Antibody.



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.