




# CDH1 Antibody

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA795479  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | P12830  |
| <b>Immunogen</b>           | Peptide sequence around aa.864~868(E-W-G-N-R) derived from Human E-Cadherin.  |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Human,Mouse,Rat   |
| <b>Specificity</b>         | The antibody detects endogenous level of total E-Cadherin protein.  |
| <b>Tested Applications</b> | ELISA,WB;WB:1:500-1:1000  |
| <b>Form</b>                | Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.                  |
| <b>Purification Method</b> | Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide. |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | CDH1; CDHE; UVO   |
| <b>Product Type</b>        | Polyclonal Antibody   |
| <b>Immunogen Species</b>   | Homo sapiens (Human)  |
| <b>Target Names</b>        | CDH1  |

## Image



Western blot analysis of extract from MCF-7 and HT29 cells using E-Cadherin(CD324) Antibody