




PECAM1 Antibody

Product Code	CSB-PA616412
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16284
Immunogen	Peptide sequence around aa.730~734(R-T-E-G-S) derived from Human CD31(PECAM1).
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total CD31 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	CD31; PECAM-1; FLJ34100; FLJ58394; PECAM1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PECAM1

Image



Western blot analysis of extract from hematopoietic stem cells using CD31(PECAM1) Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.