




FABP4 Antibody

Product Code	CSB-PA114281
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15090
Immunogen	Peptide sequence around aa.123~127(T-S-T-R-V)derived from Human FABP4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total FABP4 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	aP2; ALBP; AFABP; A-FABP;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FABP4

Image



Western blot analysis of extract from 3T3, HeLa cells using FABP4 Antibody