




FOXP4 Antibody

| | |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-PA825220 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8IVH2 |
| Immunogen | Peptide sequence around aa.7~11(S-E-T-I-R) derived from Human FOXP4. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Specificity | The antibody detects endogenous level of total FOXP4 protein. |
| Tested Applications | ELISA,WB;WB:1:500-1:1000 |
| Form | Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Purification Method | Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide. |
| Clonality | Polyclonal |
| Alias | hFKHLA; FLJ40908; FLJ44184 |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | FOXP4 |

Image



Western blot analysis of extract from JK and HeLa cells using FOXP4 Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.