




GATA4 Antibody

Product Code	CSB-PA084937
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43694
Immunogen	Peptide sequence around aa.434~438(D-S-H-G-D) derived from Human GATA-4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total GATA-4 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	MGC126629
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GATA4

Image



Western blot analysis of extract from HeLa and 293 cells using GATA-4 Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.