




GAPDH Antibody

Product Code	CSB-PA115482
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04406
Immunogen	Peptide sequence around aa.252~256(P-A-K-Y-D) derived from Human GAPDH.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total GAPDH protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	G3PD; GAPD; MGC88685
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GAPDH

Image



Western blot analysis of extract from HeLa, 293 and 3T3 cells using GAPDH Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.