




DLK1 Antibody

Product Code	CSB-PA824488
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80370
Immunogen	Peptide sequence around aa.375~379(S-K-E-A-G) derived from Human DLK1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total DLK1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	FA1; ZOG; pG2; PREF1; Pref-1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DLK1

Image



Western blot analysis of extract from HepG2 and PC12 cells using DLK1 Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.