




FGFR4 Antibody

Product Code	CSB-PA556085
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22455
Immunogen	Peptide sequence around aa.45~49(L-T-V-A-L) derived from Human FGFR4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total FGFR4 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	TKF; JTK2; CD334; MGC20292;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FGFR4

Image



Western blot analysis of extract from Rat liver and Rat brain tissue using FGFR4 Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.