




GAD2 Antibody

Product Code	CSB-PA189687
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05329
Immunogen	Peptide sequence around aa.148~152 (E-L-A-D-Q) derived from Human GAD65 (GAD2).
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total GAD65 (GAD2)protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	MGC161605; MGC161607;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GAD2

Image



Western blot analysis of extract from Rat brain, Mouse brain tissue and HeLa, 3T3 cells using GAD65(GAD2) Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.