




ITGA4 Antibody

Product Code	CSB-PA977888
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13612
Immunogen	Peptide sequence around aa.1024~1028(Y-I-N-S-K) derived from Human LPAM-1(Integrin a4, CD49d).
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total LPAM-1(Integrin a4, CD49d) protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	IA4; MGC90518; ITGA4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ITGA4

Image



Western blot analysis of extract from JK and Raji cells using LPAM-1(Integrin a4, CD49d) Antibody