




FOXJ1 Antibody

Product Code	CSB-PA979139
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92949
Immunogen	Peptide sequence around aa.413~417(D-W-A-S-V) derived from Human Foxj1 (HFH4).
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total Foxj1 (HFH4) protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	HFH4; HFH-4; FKHL13
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FOXJ1

Image



Western blot analysis of extract from SW626, HepG2, Hela and MCF-7 cells using Foxj1(HFH4) Antibody