




HMGA2 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA123352 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P52926 |
| Immunogen | Peptide sequence around aa. 6~10(E-G-A-G-Q) derived from Human HMGA2. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Specificity | The antibody detects endogenous level of total HMGA2 protein. |
| Tested Applications | ELISA,WB;WB:1:500-1:1000 |
| Form | Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Purification Method | Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide. |
| Clonality | Polyclonal |
| Alias | HMGIC; HMGI-C; STQTL9; BABL; LIPO |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | HMGA2 |

Image



Western blot analysis of extract from HeLa and 3T3 cells using HMGA2 Antibody