




RING1 Antibody

Product Code	CSB-PA205793
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06587
Immunogen	Peptide sequence around aa.35~39(I-A-V-S-P) derived from Human Ring1A.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total Ring1A protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	RNF1; RING1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RING1

Image



Western blot analysis of extract from HeLa cells using Ring1A Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.