




Postn Antibody

Product Code	CSB-PA202817
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q62009
Immunogen	Peptide sequence around aa. 125~129(R-E-E-I-E) derived from Mouse Periostin.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total Periostin protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	POSTN; OSF2; OSF-2; PN; osteoblast specific factor
Product Type	Polyclonal Antibody
Immunogen Species	Mus musculus (Mouse)
Target Names	Postn

Image



Western blot analysis of extract from Rat brain and Mouse brain Tissue using Periostin Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.