




MBP Antibody

Product Code	CSB-PA261629
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02686
Immunogen	Peptide sequence around aa.291~295(G-G-R-D-S
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total MBP(myelin basic protein) protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	MGC99675
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MBP

Image



Western blot analysis of extract from Rat brain, Mouse brain and Rat hippocampus Tissue using MBP Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.