




FOXP1 Antibody

Product Code	CSB-PA947591
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H334
Immunogen	Peptide sequence around aa. 669~673 (D-E-P-V-N)
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total FOXP1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	QRF1; 12CC4; hFKH1B; HSPC215;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FOXP1

Image



Western blot analysis of extract from JK, Raji, HeLa and 293 cells using FOXP1 Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.