




Grin2b Antibody

Product Code	CSB-PA191146
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00960
Immunogen	Peptide sequence around aa.1250-1254(N-L-Y-D-I) derived from Human NMDANR2B Subunit.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total NMDANR2B Subunit protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	NMDAR2B; NR2B;
Product Type	Polyclonal Antibody
Immunogen Species	Rattus norvegicus (Rat)
Target Names	Grin2b

Image



Western blot analysis of extract from T98G cells using NMDANR2B Subunit Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.