




FOXA2 Antibody

Product Code	CSB-PA146585
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y261
Immunogen	Peptide sequence around aa.435~439(A-S-P-L-A) derived from Human FoxA2(HNF3b).
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total FoxA2(HNF3b) protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	HNF3B; TCF3B;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FOXA2

Image



Western blot analysis of extracts from 293, HT29, HepG2 cells using FoxA2(HNF3b) Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.