




CEBPB Antibody

Product Code	CSB-PA828248
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17676
Immunogen	Peptide sequence around aa.336~340(E-P-L-L-A) derived from Human C/EBPb.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total C/EBPb protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	IL6DBP; CRP2; TCF5; LAP; NF-IL6
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CEBPB

Image



Western blot analysis of extracts from Rat brain, Mouse brain tissue and 293 cells using C/EBPb Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.