




ANXA6 Antibody

Product Code	CSB-PA153225
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08133
Immunogen	Peptide sequence around aa.6~10(Q-G-A-K-Y) derived from Human Annexin VI.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total Annexin VI protein.
Tested Applications	ELISA, WB; WB: 1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	ANX6; CBP68;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ANXA6

Image



Western blot analysis of extracts from HepG2 and HeLa cells using Annexin VI Antibody.