




JUP Antibody

Product Code	CSB-PA909861
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14923
Immunogen	Peptide sequence around aa.726~730(I-D-T-Y-S) derived from Human g-Catenin
Raised In	Rabbit
Species Reactivity	Human,Rat
Specificity	The antibody detects endogenous levels of total g-Catenin protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Clonality	Polyclonal
Alias	DP3; PDGB; PKGB; CTNNG; ARVD12
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	JUP

Image



Western blot analysis of extracts from HT-29, MCF-7 and C2C12 cells using g-Catenin Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.